ACT-VETube

ad us. vet.

In vitro Diagnostikum

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ad us. vet.

 In vitro diagnosticum

 Screening test for evaluation of
secondary hemostasis disorders of the dog and cat

INSTRUCTIONS FOR USE

1. INFORMATION ON THE TEST-KIT

TEST-KIT COMPONENTS

1 test-kit **ACT-*VE*Tube** contains:

- 10 **ACT-*VE*Tube** tubes coated with diatomaceous earth
- 1 instructions for use

STABILITY AND STORAGE

Store at
15–25°CExpiry date
– see label

APPLICATION

For veterinary use only

Lot number

In vitro diagnosticum

Do not use test-kit components from different kits, lot numbers or beyond stated expiry date.

Follow instructions for use precisely

LIABILITY

The entire risk due to the performance of this product is assumed by the purchaser. The manufacturer shall not be liable for indirect, special or consequential damages of any kind resulting from the use of this product.

2. INTRODUCTION

The Activated Clotting Time (ACT) is, like the Activated Partial Thromboplastin time (aPTT), an important coagulation marker for screening of the individual status of haemostasis, especially of the intrinsic and common path of the coagulation cascade.

Haemostasis is a complex physiological reaction to all kinds of bleedings. Bleeding disorders can be divided into primary and secondary defects or causes, respectively.

With ACT, suspicion of secondary haemostasis disorder can be confirmed quickly and easily on-site. Prolonged ACT occurs during internal bleedings in body cavities, by genetically inherited factor deficiencies like haemophilia A and B as well as by acquired causes like rodenticide intoxication, liver diseases, administration of drugs like coumarin or heparin and lack of vitamin K.

Being fast, simple and reliable, **ACT-*VE*Tube** enables the veterinarian to measure the time (in seconds) needed by a blood sample mixed with diatomaceous earth at 37°C to show first signs of blood clots.

3. INFORMATION ON THE SPECIMEN MATERIAL PATIENT PREPARATION

Usual preparations of the patient for blood withdrawal. First choice of venipuncture site (optimally with 20-gauge needle) is the Vena jugularis, alternatively also cephalic and saphenous venipuncture are possible. It is essential that a single puncture and a steady and fast flow of blood is guaranteed. Thereby, the first drops/stream must be discarded to avoid contamination of the sample with tissue factors like factor III, thromboplastin etc., because these may hasten the clotting process.

4. TEST PREPARATION

- Warm up the **ACT-*VE*Tube** to ca. 37°C (fig.1).
- Blood withdrawal from the V. jugularis (fig.2), V. cephalica antebrachii or V. saphena lateralis. Discard the first drops.

5. TEST PROCEDURE

- Add 2.0 ml native blood (freshly drawn, without anticoagulant!) from V. jugularis, V. cephalica or V. saphena into the hand-hot **ACT-*VE*Tube** (fig.3).
- Start stop-watch immediately (fig.4).
- Close **ACT-*VE*Tube** and mix by careful swinging for one time (fig.5).
- Swing the tube from horizontal to vertical every 10 seconds and watch for potential coagulation (fig.6).

6. READING OF THE TEST RESULT

When **first blood clots are visible**, stop the stop-watch and read coagulation time (ACT, fig.7).


fig.1

fig.2

fig.3

fig.4

fig.5

fig.6

fig.7

7. PRECAUTIONS FOR USERS

- Use a new **ACT-*VE*Tube** for each patient.
- The sample material must be seen as potentially infectious and disposed of accordingly, together with the used test-kit components.

NOTE:

Performing **ACT-*VE*Tube** in cats can be difficult. To obtain reliable test results, the **ACT-*VE*Tube** should be directly filled with fresh blood from a vacutainer system, preferable under sedation! Cats could be much easier tested using the SCA 2000 coagulometer for the detection of PT and aPTT values!

8. TEST PRINCIPLE

The Activated Clotting Time (ACT) is defined as the time in seconds for whole blood mixed with diatomaceous earth until the appearance of the first unmistakable clots at 37°C. Fresh native blood from V. jugularis is directly added into the **ACT-*VE*Tube** coated with diatomaceous earth, closed and mixed carefully. Activating the intrinsic cascade, diatomaceous earth will coagulate native blood in a definite time due to the coagulation status of the animal. Once first clot signs form, the test is complete and the required time (= ACT time) has to be noted.

9. INFORMATION FOR THE INTERPRETATION

ACT standard time

- Dog: 60–110 seconds
- Cat: 50–75 seconds

Prolonged ACT time

- DIC (Disseminated Intravascular Coagulation)
- advanced rodenticide or fungal poisoning, snake bite
- severe thrombocytopenia (< 10.000/μl blood)
- diverse thrombocytopathies
- inhibitors preventing clot formation (heparin, coumarin, uraemia, a.s.o.)
- severe factor deficiencies (sensitivity < than aPTT!) – intrinsic or common (XII, XI, IX, VIII, X, V, II, I)

Measuring the time between the first clot signs and formation of a solid clot involving the whole sample could give additional information. This period is often increased in dogs with coagulopathy, in dogs with normal coagulation it is short (seconds).