# FASTest® BCV Strip

ad us.vet.

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# FASTest® BCV Strip

ad us.vet.

Test-kit for the detection of  
Bovine Coronavirus  
antigen in calf feces

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### 3. Test-Kit Components

#### 1 Test-kit FASTest® BCV Strip:

- 2 or 10 Dipstick Tests coated with monoclonal antibodies.
- 2 or 10 Sample Tubes with 3,0 ml buffer diluent.
- 1 Instructions for use.

### 4. Storage and Shelf Life

- Store at room temperature (+15°C - + 25°C).
- Stored correctly the product can be kept up to the expiry date.
- Avoid the Test kit being subjected to excess heat or freezing.

### 1. SPECIMEN COLLECTION AND DILUTION

- Open the sample tube and take out the sample spoon, provided with the kit.
- Take a little portion (0,4g) of feces with the spoon and put it back into the sample tube. Screw on the cap tightly avoiding leakage. (fig.1)
- If possible let the test material dissolve in the sample tube.
- Shake the tube until the sample has been dissolved into the buffer diluent. (fig.2)
- Leave the test tube until the large particles have settled to the bottom of the tube.

### 1. Introduction

Bovine Coronavirus infection is one of the most frequent cause of Gastroenteritis by calves and related species. Coronavirus infections in calves is one of the most common germ of neonatal diarrhea. The viruses breeds in the small intestine and are secreted with feces. After secretion the bovine Coronavirus remains for a long time in the environment.

Coronavirus infections show fluid and mucus diarrhea, mostly over a short period of time. In adult animals, Coronavirus infections have a sub-clinical course. A passive transfer of Coronavirus maternal antibodies have been demonstrated. Stress situations and secondary viral or bacterial infections concurrent to bovine Coronavirus infections could prove the fatal to the animal.

### 5. Information on the Test Sample Material

Only feces can be examined using the **FASTest® BCV Strip**. The test sample material can be stored for up to 14 days at 2°-8°C. If a longer storage period is required then the material must be prepared as follows:

Use a pipette to extract a sufficient quantity of the test fluid from the surface of the tube and store this in a suitable container ( e.g. microcentrifuge tube ) at -25°C.

### 6. Special Information

- FOR VETERINARY USE ONLY!
- The test can be used with calf or bovine feces.
- Do not use the test once the expiry date has passed.
- Do not use reagents from different kits.
- The sample material must be seen as potentially infectious.
- Follow instructions for use precisely.

### 2. TEST PROCEDURE

- Remove the dipstick from its pouch by tearing along the split.
- Label the dipstick with a patient name or ID number.
- Introduce the dipstick vertically with the arrows pointing downwards into the sample tube with the feces-buffer-fluid. (fig.3)
- Remove the dipstick from the sample tube with the feces-buffer-fluid after 1 Minute or longer. (fig.4)
- Place the dipstick on a flat surface.

### 3. INTERPRETATION OF TEST RESULTS

Read the results 5 minutes after the dipstick has been placed in the sample tube with feces-buffer-fluid. The test result will be stable over a long term.

### 2. Test Principle

The **FASTest® BCV Strip** is based on a rapid immunochromatographic test principle. The test membrane contains two unique monoclonal antibodies. One of these monoclonal antibodies is bonded to colloidal gold particles, the second is immobilized in the test zone of the strip. If the stool sample extract contains bovine Coronavirus antigen these will form an antigen-antibody complex with the goldparticles. As this complex migrates along the TEST zone of the membrane to the immobilized capture antibody, a pink-purple line is formed indicating a positive test result.

The remaining conjugate migrates to a second antibody on the CONTROL zone of the test strip forming a pink-purple line. This indicates proper performance of the test.

### LIABILITY

**The entire risk due to the performance of this product is assumed by the purchaser. The manufacturer shall not be liable for indirect, special or consequential damages of any kind resulting from use of this product.**

### POSITIVE TESTRESULT

In addition to the pink line in the CONTROL zone, a clearly visible pink-purple line also appears in the TEST zone. The feces sample contains bov. Coronavirus. (fig.5)

### NEGATIVE TESTRESULT

Only the pink-purple line appears in the CONTROL zone. No visible line is seen in the TEST zone. The feces sample contains no bov. Coronavirus. (fig.6)

### INCONCLUSIVE TESTRESULT

If there are no pink-purple lines, neither in the TEST zone nor the CONTROL zone of the dipstick, visible after 10 minutes. The test should be repeated using a new dipstick.

### POSITIVE TESTRESULT

### NEGATIVE TESTRESULT