# FASTest® CDV Ab

ad us. vet.

In vitro Diagnostikum

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ad us. vet.

*In vitro* diagnosticum

Test-kit for the semiquantitative detection of  
Canine Distempervirus antibodies  
in whole blood, plasma or serum of the dog

## INSTRUCTIONS FOR USE

## 1. INFORMATION ON THE TEST-KIT

### TEST-KIT COMPONENTS

1 test-kit **FASTest® CDV Ab** contains:

- 2 or 10 test cassettes coated with recombinant antigens
- 2 or 10 buffer diluent tubes **A** with 0.25 ml buffer diluent each
- 2 or 10 disposable inoculation loops
- 2 or 10 disposable plastic pipettes
- 1 instructions for use

### STABILITY AND STORAGE

Store at  
15–25°C  
15–25°C

Expiry date  
– see label

### APPLICATION

For veterinary use only

Lot number

*In vitro* diagnosticum

Do not use test-kit components from different kits, lot numbers or beyond stated expiry date.

Follow instructions for use precisely

### LIABILITY

The entire risk due to the performance of this product is assumed by the purchaser. The manufacturer shall not be liable for indirect, special or consequential damages of any kind resulting from the use of this product.

## 2. INTRODUCTION

Antibodies are basic modules of the humoral immune response. They are passed by passively via the colostrum as so-called maternal antibodies (mAb) onto the yet immunocompetent newborns or induced actively by natural field infection or vaccination. The antibody titre is varying individually in each animal, depending on multiple factors. The titre can persist over an extended period of time, partially lifelong, in efficient protection concentration (= reliable immunity by protective antibodies) or can fall below the efficient protection concentration (non-reliable immunity) in the course of time.

Depending on the level of individual antibody titre, the veterinarian is able to decide fast and reliable the necessity of vaccination or non-vaccination due to following questions:

Individual vaccination point of the breeding bitch: In problematic breedings, the determination of antibody status of the female makes sense during pregnancy to decide whether a booster vaccination before birth is necessary or to find the optimal primary vaccination time of the puppies.

There is a critical stage (so-called immunity gap) in puppies, especially in the first 12 weeks. During this stage the concentration of mAb could be high enough to inactivate the vaccinating virus but also too low to protect from field infection. Therefore it is important to find the individual primary vaccination point for each puppy to guarantee an appropriate protection.

For the determination of antibody status of the whole litter, it is possible to determine the antibody status of only one puppy, representative for the other puppies (so-called "fraternal antibody titre"). Here, the balanced colostrum assumption or development of all puppies is absolutely necessary.

By determination of the actual antibody status, an individual decision of the necessity of booster vaccination of the puppy or the adult animal can be made.

Being fast, safe and reliable, for pet owner and breeder these important questions can be answered practically by **FASTest® CDV Ab**. This enables the veterinarian an appropriate and customized vaccination diagnostics and strategy, adapted to dog and pet owner.

## 3. INFORMATION ON THE SPECIMEN MATERIAL

Approximately 1 µl (1 inoculation loop) 15–25°C warm whole blood (WB, native blood with anticoagulant), plasma (P) or serum (S) are needed. Native blood without anticoagulant should not be used due to potential micro agglutination (e.g. migration delay on the membrane, unspecific reaction)!

Mix the sample material well before use!

Non-cooled (15–25°C), WB, P, S and E should be tested within 4 hours! At 2–8°C, WB, P and S can be stored up to 4 days. Serum and/or plasma samples can be permanently stored at minimum –20°C.

Keep in mind that the sample material, as well as all used test-kit components, should have reached room temperature at the time of application.


Endogeneous and exogeneous interfering substances of the sample (e.g. albumin, fibrinogen, lipids, CRP, heterophilic antibodies, especially type IgA, as well as viscosity, pH-value and excess EDTA) can cause interferences (matrix effects) that can influence the target measurement. These can lead to an impaired LF and/or unspecific reactions on T and C1 or C2.

## 4. SPECIMEN COLLECTION AND PREPARATION

- a. Pick up 1 µl of specimen using the disposable inoculation loop by dipping the circular end of the loop into the specimen (fig.1).
- b. Open the yellow cap of the buffer diluent tube **A** and add the amount of specimen (1 µl) into the buffer diluent of the buffer diluent tube **A** (fig.2).
- c. Gently stir the loop with specimen into the buffer diluent to ensure adequate mixing of the specimen into the buffer diluent. (fig.2). The remove the loop and close again the buffer diluent tube **A**. Discard the loop.

## 5. TEST PROCEDURE

1. Remove the test cassette from its foil pouch shortly before use. Place it on a flat surface.
2. Open the buffer diluent tube **A** containing the homogeneously mixed buffer-sample mixture (BSM). Take the disposable plastic pipette and suck off the whole content and drop the whole BSM into the round sample window of the test cassette (fig.3).
3. Add 1 additional drop of buffer diluent into the round sample window if there is no beginning LF visible within 1 minute after adding the buffer diluent.

## 6. READING OF THE TEST RESULT

Read the test result 20 minutes after addition of the BSM into the round sample window.

### POSITIVE TEST RESULT

(prolonged vaccination interval possible)

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
### NEGATIVE TEST RESULT

(Distempervirus vaccination recommended)

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### INVALID TEST RESULT

No **T** line and no **C** lines are visible. The test should be repeated using a new test cassette.

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## 7. PRECAUTIONS FOR USERS

- Label sample material and associated test cassette to ensure a precise assignment.
- Use a new inoculation loop, buffer diluent tube and pipette for each sample.
- The buffer diluent contains low concentrations of toxic sodium azide as a preservative, therefore avoid skin contact and/or ingestion.
- The sample material must be seen as potentially infectious and disposed of accordingly, together with the used test-kit components.

## 8. TEST PRINCIPLE

The **FASTest® CDV Ab** is based on an immunochromatographic "sandwich principle".

In the sample pad are recombinant CDV antigens, in the conjugate pad two different mobile mono or polyclonal antibodies bound to gold particles. The membrane shows three lines: the CONTROL lines **C1** and **C2** and the TEST line **T**. Each line is coated with polyclonal antibodies.

The anti-CDV antibodies of the sample first react with the recombinant CDV antigens of the sample pad, second with the mobile monoclonal gold labeled antibodies of the conjugate pad. During the following "lateral flow" (LF) along the nitrocellulose membrane, these antigen-antibody complexes are captured by fixed polyclonal antibodies forming a pink-purple test line **T**. The colour intensity of **T** can vary depending on the anti-CDV antibody concentration of the sample.

In parallel, also the mobile polyclonal gold labeled antibodies of the conjugate pad migrate along the membrane (LF) and react with the fixed polyclonal antibodies of **C1** and **C2** forming pink-purple coloured CONTROL lines. The formation of **C1** and **C2** indicates a correct test procedure.

Evaluation of **FASTest® CDV Ab** is done by comparison of the colour intensities of **T** with **C1** and **C2**.

The threshold titre (sustainable immunity or not) of **FASTest® CDV Ab** (1:80) is adjusted by Golden Standard Test (hemagglutination inhibition test).

## 9. INFORMATION FOR THE INTERPRETATION

- The interpretation of the test result should always be based on anamnestic and clinical data as well as the therapy and prophylaxis possibilities.
- Any non-described colour or contour variation of **TL** and **CL** (e.g. greyish, shadow-like lines) has to be considered as unspecific reactions and therefore as negative test result.
- Due to anticoagulated whole blood and/or red hemoglobin background of the test membrane, caused by hemolytic blood samples, the visibility of **T**, especially in case of weak positive samples, could be from worse to not visible.
- Any coloured lines appearing after 30 minutes do not have any diagnostic value.
- The **FASTest® CDV Ab** only detects the presence or absence of anti-CDV IgG antibodies in the specimen and should not be used as the sole criterion for the diagnosis of CPV infection in dogs.

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