FASTest® GIARDIA Strip

ad us. vet.

In vitro Diagnostikum

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ad us. vet.

In vitro diagnosticum

Test-kit for the qualitative detection of *Giardia duodenalis* surface antigens in feces of the dog and cat

INSTRUCTIONS FOR USE

The test is approved for dogs and cats only.

1. INFORMATION ON THE TEST-KIT

TEST-KIT COMPONENTS

1 test-kit **FASTest® GIARDIA** Strip contains:

- 2 or 10 dipsticks coated with monoclonal antibodies
- 2 or 10 sample tubes with 2.0 ml buffer diluent each
- 1 instructions for use

STABILITY AND STORAGE

Store at
15–25°C
15–25°C

Expiry date
– see label

APPLICATION

For veterinary use only

LOT

Lot number

In vitro diagnosticum

Do not use test-kit components from different kits, lot numbers or beyond stated expiry date.

Follow instructions for use precisely

TL – TEST line, **CL** – CONTROL line, **LF** – Lateral flow

LIABILITY

The entire risk due to the performance of this product is assumed by the purchaser. The manufacturer shall not be liable for indirect, special or consequential damages of any kind resulting from the use of this product.

2. INTRODUCTION

Giardia is known to be one of the most common enteric parasites in pocket pets, pets, farm and wild animals as well as in humans (zoonosis) world-wide.

G. duodenalis occurs in varying genotypes (genotypes A to G), which differ in their infection and their host spectrum. Types A and B have zoonotic potential.

Newborns and young animals are mostly affected. Prevalences vary in cats and dogs, depending on age (>70% under 1 year), husbandry (10% in single husbandry up to 100% in breedings and animal shelters) and immune status.

Transmission (direct contact, by contaminated food, water, objects, grooming and vectors like flies etc.) happens fecal-orally by ingestion of highly infectious and very resistant cysts being discharged by other animals or humans. Only five to seven cysts are enough to cause an infection.

G. duodenalis has an asexual life cycle. In the duodenum of the infected animals, two so-called trophozoites, emerge from the incorporated cysts (excystment). These multiply by duplication and attach via suckers to the duodenal surface. Free trophozoites turn into their last forms, the cysts (encystment), especially in the ileum. These are excreted in large amounts (10⁷/g feces) and mostly intermittent, i.e. not with every defecation. The prepatent period averages 5 to 16 days.

The main symptom of giardiasis is diarrhoea, more or less intensive, that can run from symptomatic (acute, chronic, self-limiting, periodic-intermittent or continuous) to asymptomatic. Independent on the progression, cysts and/or trophozoites can be egested (primarily with strong diarrhoea).

Giardia cysts can be differentiated from cysts of different coccidia species only by microscopical experienced people. This is similarly true for *Giardia*- and *Trichomonas foetus* trophozoites. For that reason, the modern aetiological coproscopy via **FASTest® GIARDIA** Strip should be preferred compared to the microscopical proof.

For epidemiological reasons, all animals, clinical symptomatic and clinical asymptomatic, have to be tested with **FASTest® GIARDIA** Strip. This enables the veterinarian in the clinic to introduce a specific treatment as well as a broad prophylaxis.

3. INFORMATION ON THE SPECIMEN MATERIAL

Due to the normally inhomogeneous or nest-like dissemination of antigens in the feces, the specimen material has to be mixed up homogeneously (spatula, vortex-mixer) before sampling.

For the test, the required amount of feces as described in issue 4b/Specimen collection and preparation, is needed. The amount depends on the consistency of the sample. Use the attached spoon.


Non-cooled (15–25°C), the sample should be tested within 4 hours! At 2–8°C, the sample can be stored up to 4 days, permanently at minimum –20°C.

Keep in mind that the sample material, as well as all used test-kit components, should have reached room temperature at the time of application.

Endogeneous and exogeneous interfering substances of the sample (e.g. proteases, mucosa components, blood, but also viscosity, pH-value as well as grass and cat litter) can cause interferences (matrix effects) that can influence the target measurement. These can lead to an impaired LF and/or unspecific reactions on the TL and CL.

4. SPECIMEN COLLECTION AND PREPARATION

- Open the sample tube with the buffer diluent.
- Mix the feces sample homogeneously (applicator, vortexer). Then mix the required sample volume (fig.1):



pulpy-compact: 1, fluid-watery 2 to max. 3 level spoons of feces steadily into the buffer diluent.

- Close sample tube tightly and rotate it easily to get the mixture as homogeneous as possible (fig.2).
- For sedimentation of gross feces particles place the sample tube on a flat and horizontal surface for 1–5 minutes.

5. TEST PROCEDURE

- Remove the dipstick from its foil pouch shortly before use.
- Introduce the dipstick vertically and with the arrows pointing downwards into the sample tube for at least 1 minute. The liquid level must not exceed the red plastic cover with the white arrows (fig.3).
- Remove the dipstick from sample tube soonest the sample-buffer mixture (SBM) has reached the CL. If so, the blue CL will appear slowly but surely (fig.4). If the CL will not appear after 5–10 minutes, a new SBM must be prepared and sedimented for at least 5 minutes. The dipstick must be held only in the supernatant until the LF has reached the CL.
- Place the dipstick on a flat and horizontal surface (fig.4).

6. READING OF THE TEST RESULT

Read the test result after **5 minutes**. Positive test results may be observed earlier, depending on the concentration of antigen in the sample.

POSITIVE TEST RESULT (fig.5)

A weak to strong and well-defined **red TEST line** and a **well-defined blue CONTROL line** appear.

NEGATIVE TEST RESULT (fig.6)

Only a **well-defined blue CONTROL line** appears. This line indicates, irrespective of its intensity, that the test has been performed properly.

INVALID TEST RESULT

Only a weak to strong and well-defined red TEST line or no line at all appears. The test should be repeated using a new dipstick.

fig.5

POSITIVE TEST RESULT

fig.6

NEGATIVE TEST RESULT

7. PRECAUTIONS FOR USERS

- Label sample material and associated sample tube to ensure a precise assignment.
- Use a new sample tube for each sample.
- The buffer diluent contains low concentrations of toxic sodium azide as a preservative, therefore avoid skin contact and/or ingestion.
- The sample material must be seen as potentially infectious and disposed of accordingly, together with the used test-kit components.

8. TEST PRINCIPLE

The **FASTest® GIARDIA** Strip is based on latest rapid immunochromatographic technique. Surface antigens of intact or fractionated *Giardia duodenalis* cysts and/or trophozoites will react at the conjugate pad with mobile antibodies bound to red latex particles. Migrating ("lateral flow", **LF**) along the nitrocellulose membrane, these specific antigen-antibody complexes are bound by fixed monoclonal anti-*Giardia* antibodies (mAbs) producing a more or less red **TEST line (TL)**. These mAbs guarantee a high level of specificity for the aetiological detection of *G. duodenalis* antigens. A correct test procedure will be indicated by a second, intense blue **CONTROL line (CL)**. The intensity or width of the TL depends on the concentration of *Giardia* antigens in the introduced amount of sample. The **FASTest® GIARDIA** Strip does not rely on the presence of intact cysts and/or trophozoites.

9. INFORMATION FOR THE INTERPRETATION

- The interpretation of the test result should always be based on anamnestic and clinical data as well as the therapy and prophylaxis possibilities.
- Any non-described colour or contour variation of TL and CL within the indicated incubation time or after more than 10 minutes (e.g. greyish, shadow-like lines) has to be considered as unspecific reaction and therefore as negative test result.
- TL can vary both in intensity (from weak to intense red) and width. Therefore, any red line appearing within the required incubation time is to be interpreted as a positive test result.
- A negative test result does not exclude a persistent *Giardia duodenalis* infection due to the intermittent shedding (not with every defecation!) of cysts or trophozoites. In case of an ongoing diarrhoea and a negative test result, the test should be repeated after one week using a serial feces sample (one feces every day for three days or three feces samples consecutively). Only a multiple negative feces test rules out a *Giardia* infection. Any red test line, despite of intensity and/or width is to be interpreted as a positive test result.
- Drug therapy can lead initially to an increase of cyst and/or trophozoite shedding. Reinfections due to the prepatence could also appear. Both could result in positive test results. First of all, therapy control should be based on clinical symptoms (decrease of diarrhoea in quantity and quality), not only of a single test result.
- Be aware, every dog or cat tested positive is considered potentially infectious for animals and humans (zoonosis!), especially for kids!