FASTest® IgG bovine

ad us. vet.

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

ad us. vet.

In vitro diagnosticum

Test-kit for the semiquantitative detection of IgG antibodies in native blood, whole blood, plasma and serum of cattle

INSTRUCTIONS FOR USE

A-6912 Hörbranz – AUSTRIA

1. INFORMATION ON THE TEST-KIT

TEST-KIT COMPONENTS

- 1 Test-kit **FASTest® IgG bovine** contains:
- 2 or 10 test cassettes coated with polyclonal antibodies
 - 1 dropper bottle **A** with 1.0 ml or 3.0 ml buffer diluent
 - 4 or 20 disposable plastic pipettes
 - 2 or 10 disposable sample tubes
 - 1 instructions for use

EXPLANATION OF THE SYMBOLS

STORAGE AND SHELF LIFE

- Store at room temperature (15–25 °C).
- Stored correctly, the product can be used up to the expiry date.
- Avoid the test-kit being subjected to excessive heat or freezing.

LIABILITY

The entire risk due to the performance of this product is assumed by the purchaser. The manufacturer shall not be liable for indirect, special or consequential damages of any kind resulting from the use of this product.

2. INTRODUCTION

An effective immunity is based on optimal levels of immunoglobulins (IgG, especially IgG1). This guarantees an optimal breeding success, especially in the periparturial stage and in new born calves, and is of economic importance to the livestock industry.

Due to the special placenta conditions, newborn calves often show no to hardly noteworthy IgG. Therefore the most important base for an immunoprophylaxis is an adequate supplementation with IgG.

IgG transfer mainly is based on the ingestion and resorption of colostrum (adequate passive transfer) directly after birth up to maximal 24 h later. Failure of passive transfer (FPT) may be caused by inadequate suckling (reduced vitality, neonatal respiratory depression a.s.o.), way of colostrum application, very low levels of IgG in the colostrum, inadequate absorption of IgG or environmental stress. The consequences of an immune deficiency syndrome (IDS) are weak calves with a higher vulnerability for infectious newborn diseases like enzootic pneumonia of calves, diarrhoea (neonatal calf diarrhoea; calf scours) and other septicaemic diseases.

The incidence for cows (dystocia, metritis, mastitis a.s.o.) in the periparturial period is as high as in no other farm animal. Latest investigations show following IgG1 concentrations in healthy cows with normal parturition and physiologic postparturial period:

- Ø 38 mg/ml at 7th month of gestation
 - Ø 15 mg/ml physiologically low concentration periparturial
 - Ø 30 mg/ml increasing from 4th month post partum
- If the intraparturial IgG concentration is already significantly decreased (< 15 mg/ml), an increased incidence of genital associated diseases will be seen in the periparturial period.

The optimal test time slot using **FASTest® IgG bovine** is between 24 to 48 h (max up to 7 days) post natum in calves and from 3rd to 7th day post partum in cows. **FASTest® IgG bovine** enables the veterinarian on-farm and without any technical costs to confirm (**FASTest® IgG bovine**: IgG ≤ 12 mg / ml) or to exclude (**FASTest® IgG bovine**: IgG > 12 mg / ml) a suspicion for FPT/IDS in the cow as well in the calf.

3. INFORMATION ON THE SPECIMEN MATERIAL

For the test **30-35 µl (1 drop using the disposable pipette) of freshly obtained and 15–25 °C warm native blood** (NB, without anticoagulant; should be tested immediately after blood withdrawal!) or **whole blood** (WB, with anticoagulant), **plasma** (P) or **serum** (S) is required.


Non-cooled (15–25 °C) NB should be used immediately, WB P and S latest within 4 hours after blood withdrawal! At **2–8 °C**, WB, P and S can be stored up to 4 days. **Serum- and/or plasma samples** can be permanently stored at **minimum -20 °C**.

The ratio anticoagulant to necessary blood volume must be ensured to avoid migration delay during the “lateral flow” (e.g. by microclots or unspecific reactions based on excess of anticoagulant) in the area of the T and C line.

4. SAMPLE PREPARATION (PREDILUTION)

ONLY PLASMA/SERUM!

- Take the plastic pipette and place **1 (one) drop (30-35 µl) of plasma or serum** into the associated sample tube (fig.1). Hold the pipette vertically (fig.2).
- Add **2 (two) drops (ca. 60-70 µl) buffer diluent** from the dropper bottle **A** (vertically) into the sample tube (fig.2).
- Close the sample tube and mix the sample-buffer fluid homogeneously.


5. TEST PROCEDURE

- Remove the test cassette shortly before use from its foil pouch, label it with patient name or ID number and make sure that the test cassette is on a flat surface.
- Take a new plastic pipette (hold it vertically!) and place **a NATIVE/WHOLE BLOOD b PLASMA/SERUM** 1 drop (ca. 30-35 µl) of 1 drop (ca. 30-35 µl) of the non-prediluted native homogeneously mixed sample-blood/whole blood sample-buffer fluid into the SAMPLE window “S” of the test cassette (fig.3a/b).

- After the drop has been sucked in completely, add **1 (one) drop of buffer solution (ca. 30-35 µl) of the dropper bottle A** (hold vertically!) into the SAMPLE window “S” of the test cassette (fig.4).
4. Read the test result after an incubation time of **5 minutes** in the RESULT window “C-T”.

7. PRECAUTIONS FOR USERS

- FOR VETERINARY USE ONLY!
- Do not remove any test-kit components (room temperature!) from their individually sealed pouches until shortly before their use.
- Do not use test-kit components from different kits, lot numbers or beyond the stated expiry date.
- The sample material must be seen as potentially infectious and disposed of accordingly, together with the used test-kit components.
- Follow instructions for use precisely.

8. TEST PRINCIPLE

FASTest® IgG bovine is designed on latest lateral flow technology. Based on an immunochromatographic sandwich principle, bovine IgG antibodies in the sample indirectly compete with highly specific polyclonal, mobile and gold labeled anti-bovine IgG antibodies for the binding of fixed bovine antibodies of the TEST line “T”.

Due to the competitive test principle, the intensity of the TEST line “T” is reciprocally proportional to the concentration of bovine IgG in the sample. Therefore, in case of **low IgG concentrations (high IgG concentrations) only few (a great many)** binding sites of the mobile anti-bovine IgG antibodies of the sample are blocked. These will bind in **great (in small to zero decreasing)** quantities to fixed bovine IgG antibodies at the TEST line “T” creating a **strong (weak to no)** intensive TEST line “T”. Remaining and free mobile anti-bovine IgG antibodies migrate further to the CONTROL line “C”. Depending on their concentration they create a more or less strong to very weak pink-purple CONTROL line “C”.

6. READING OF THE TEST RESULTS

Read test result **5 minutes** after dropping the prediluted (P/S) or non-diluted (NB/WB) sample and the buffer were dropped into the round sample window of the test cassette.

The evaluation of the test results is designated by reading the various colour intensities of the **TEST and CONTROL lines (fig. 5a-g)**. The competitive test principle implies: the **weaker the intensity of the TEST line, the higher the IgG concentration in the sample.**

IgG CONCENTRATION ≤ 12 mg/ml (fig.5 a, b, c, d)

The intensity range of the pink-purple TEST line “T” as well as the pink-purple CONTROL line varies from very strong to very weak.

IgG-CONCENTRATION > 12 mg/ml (fig.5 e, f, g)

The intensity range of the pink-purple TEST line “T” varies from very weak to not visible (lack of T line). The intensity range of the pink-purple CONTROL line “C” varies from very weak to extremely weak.

The absence of the C line could also appear due to extremely high concentrated IgG samples.


INVALID TEST RESULT

The test is invalid when no lateral flow (no pink-purple coloured migration of the sample-buffer fluid) appears. Reasons could be a missing predilution or using too low sample and/or buffer volume, respectively.

The interpretation of the test results should be done in the framework of anamnesis, clinical symptoms, therapy and prophylaxis measurements.

fig.5

IgG CONCENTRATION ≤ 12 mg/ml

IgG CONCENTRATION > 12 mg/ml