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FASTest® NEOSPORA**caninum**
ad us. vet.*In vitro* diagnosticumTest-kit for the qualitative detection of
anti-*Neospora caninum* antibodies in whole blood,
plasma or serum of the dog and cattle**INSTRUCTIONS FOR USE****1. INFORMATION ON THE TEST-KIT****TEST-KIT COMPONENTS**1 test-kit **FASTest® NEOSPORA** caninum contains:

- 2 or 10 test cassettes coated with recombinant *Neospora caninum* antigens
- 1 dropper bottle **A** with 1.0 ml or 3.0 ml buffer diluent
- 2 or 10 disposable plastic pipettes
- 1 instructions for use

STABILITY AND STORAGE
 Store at
15–25°C
15–25°C

 Expiry date
– see label
APPLICATION

For veterinary use only

Lot number

*In vitro* diagnosticum

Do not use test-kit components from different kits, lot numbers or beyond stated expiry date.

Follow instructions for use precisely

LIABILITY

The entire risk due to the performance of this product is assumed by the purchaser. The manufacturer shall not be liable for indirect, special or consequential damages of any kind resulting from the use of this product.

2. INTRODUCTION

Neospora caninum plays an important epidemiological role in dogs and cattle. The dog (esp. watch dogs, stray dogs) is a definitive host and excretes infectious oocysts with the feces. The intermediate host, especially cattle, but also goats, sheep and horse, gets infected via grazing land and/or water contaminated with oocysts (horizontal infection) and/or intrauterine via already infected mothers (vertical infection).

N. caninum plays an important role in abortion in cattle worldwide. Characteristics are accumulating abortions in all states of gestation, dead births and weak calves.

The dog can act as an intermediate host as well and therefore can fall sick with Neosporosis. In dogs, the symptoms are especially focused on neurological disorders: paresis/paralysis of the hind-limbs, later also of the fore-limbs, as well as polyomyelitis, radiculitis and encephalomyelitis. Also, muscular atrophy, hyperextension, hyperaesthesia and dysphagia can occur. Additionally, hepatitis, pneumonitis, myocarditis and ulcerative dermatitis can appear. In older dogs *Neospora* infection usually is asymptomatic! Puppies become clinically conspicuous at the age from 3 to 9 weeks up to one year. Early diagnostics and therefore specific therapy are essential for the prognosis. Due to recent studies, there seems to be a predisposition of male dogs to *N. caninum*.

Dog: Due to the short excretion period and the low amount of oocysts in dog feces, the detection of antibodies using **FASTest® NEOSPORA** caninum becomes very important for the diagnosis of a Neosporosis.

Cattle: Due to horizontal, but particularly to the economically more important vertical placental transmission of infection onto the offspring, with **FASTest® NEOSPORA** caninum suspicious stocks should be tested for *N. caninum* antibodies.

3. INFORMATION ON THE SPECIMEN MATERIAL

Approximately **50 µl** (2 drops of attached plastic pipette) **15–25°C warm whole blood (WB, with anticoagulant)** or **25 µl** (1 drop of attached plastic pipette) **15–25°C warm plasma (P) or serum (S)** are needed. Mix the sample material well before use!

Non-cooled (**15–25°C**), WB, P and S should be tested within 4 hours! At **2–8°C**, WB, P and S can be stored up to 4 days. **Serum and/or plasma samples** can be permanently stored at minimum **–20°C**.

Keep in mind that the sample material, as well as all used test-kit components, should have reached **room temperature** at the time of application.



Endogeneous and exogeneous interfering substances of the sample (e.g. albumin, fibrinogen, lipids, CRP, heterophilic antibodies, especially type IgA, as well as viscosity, pH-value and excess EDTA) **can cause interferences** (matrix effects) that can influence the target measurement. These can lead to an impaired LF and/or unspecific reactions on B and C.

4. SPECIMEN COLLECTION AND PREPARATION

No specimen preparation necessary.

5. TEST PROCEDURE

1. Remove the test cassette from its foil pouch shortly before use. Place it on a flat surface.
2. Add **2 drops of anticoagulated whole blood (ca. 50 µl)** or **1 drop of plasma / serum (ca. 25 µl)** with the attached plastic pipette into the sample window A of the test cassette (hold the pipette vertically, fig.1).
3. Hold the dropper bottle **A** vertically and express **4 (four) drops (160–200 µl) of buffer diluent** into the sample window A of the test cassette (fig.2).
4. Add 1 additional drop of buffer diluent into the sample window A if there is no beginning LF visible within 1 minute after adding the buffer diluent.

**6. READING OF THE TEST RESULT**

Read the test result **15 minutes** after the four drops have been added into the sample window A.

POSITIVE TEST RESULT (fig.3)

A weak to strong and **well-defined pink-purple coloured TEST line and CONTROL line** appear.

NEGATIVE TEST RESULT (fig.4)

Only a **well-defined pink-purple CONTROL line** appears. This line indicates, irrespective of its intensity, that the test has been performed properly.

INVALID TEST RESULT

Only a weak to strong and well-defined pink-purple TEST line or no line at all appears. The test should be repeated using a new test cassette.

fig.3
POSITIVE TEST RESULT

fig.4
NEGATIVE TEST RESULT

**7. PRECAUTIONS FOR USERS**

- Label sample material and associated test cassette to ensure a precise assignment.
- Use a new pipette for each sample.
- The buffer diluent contains low concentrations of toxic sodium azide as a preservative, therefore avoid skin contact and/or ingestion.
- The sample material must be seen as potentially infectious and disposed of accordingly, together with the used test-kit components.

8. TEST PRINCIPLE

The **FASTest® NEOSPORA** caninum is based on an immunochromatographic “sandwich principle” technique detecting specific anti-*Neospora caninum* antibodies in the whole blood, plasma or serum of the dog or cattle. The anti-*Neospora caninum* antibodies of the sample are migrating (“lateral flow”, **LF**) along the nitrocellulose membrane and bind to fixed recombinant *Neospora caninum* antigens conjugated with gold particles forming a pink-purple coloured TEST line (**B**). A correct test procedure will be indicated by a second pink-purple CONTROL line (**C**).

9. INFORMATION FOR THE INTERPRETATION

- The interpretation of the test result should always be based on anamnestic and clinical data as well as the therapy and prophylaxis possibilities.
- Any non-described colour or contour variation of B and C (e.g. greyish, shadow-like lines) has to be considered as unspecific reactions and therefore as negative test result.
- The proof of anti-*Neospora caninum* antibodies, together with anamnesis and clinic shows with a high likelihood that *N. caninum* can be considered as cause of the acute disease.
- In asymptomatic animals and antibody positive animals, one should assume that they could have been infected with *N. caninum* and, as potential carriers, can transmit *N. caninum* vertically on their offspring.
- Additionally, the titre amount or the titre increase (within two tests in an interval of 2–3 weeks) can be defined via indirect immunofluorescence. A proof of the pathogen can be done by PCR in foetal fluids or tissue samples of deceased whelps or in the liquor/muscle biopsies of living animals.