MYKODERMOASSAY DTM

ad us. vet.

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MYKODERMOASSAY DTM

ad us. vet.

Culture medium for the qualitative detection of veterinary relevant dermatophytes

INSTRUCTIONS FOR USE

6912 Hörbranz – AUSTRIA

1. INFORMATION ON THE TEST-KIT

TEST-KIT COMPONENTS

1 Test-kit MYKODERMOASSAY DTM contains:

- 12 sample tubes with dermatophyte culture medium
- 1 instructions for use

STABILITY AND STORAGE

Store at
2–8°C

Expiry date
– see label

APPLICATION

For veterinary use only

Lot number

In vitro diagnosticum

Do not use test-kit components from different kits, lot numbers or beyond stated expiry date.

Follow instructions for use precisely

LIABILITY

The entire risk due to the performance of this product is assumed by the purchaser. The manufacturer shall not be liable for indirect, special or consequential damages of any kind resulting from the use of this product.

2. INTRODUCTION

Dermatophytoses/ringworm belong to the most frequent infectious dermatoses in pocket pets, pets and farm animals, but also in humans (zoonosis).

They are caused by dermatophytes, filamentous fungi using keratin (skin, hair and claws) as carbon source.

The clinically most relevant species are *Microsporum* (*M. canis*, *M. gypseum*, *M. persicolor*), *Trichophyton* (*T. verrucosum*, *T. mentagrophytes*) and *Epidermophyton*. Beside age and immunosuppression, familiar, breeding (especially persian cats) and keeping conditions (breeding, animal shelter, hunting dog, multiple species keeping), travelling, lactation (transmission of infection to puppies) as well as e.g. ectoparasite based diseases and debilitated animals play an important role in developing a ringworm disease. Warm and humid climate is an additional trigger.

In case of clinical suspicion of an ongoing dermatophytosis (spotted, patchy areas of alopecia, often non-pruritic), establishment of a mycological culture using dermatophyte specific media is known to be the most reliable technique.

The MYKODERMOASSAY DTM is a classical dermatophyte test medium in sample tubes with tilted agar. It ensures a fast evaluation of the clinically suspected diagnosis through colour change of the DTM medium.

This enables the veterinarian in suspicious cases to identify a dermatophytosis and initiate a specific therapy.

3. INFORMATION ON THE SPECIMEN MATERIAL

Sampling should be done optimally before an antibiotic therapy!

Taking an optimal sample (amount/purity degree) is the most critical step for the growth of potential dermatophytes. The more hair, dandruff and/or scrabs are used, the less the likelihood to get a false positive fungal culture result.

The growth and therefore the identification of the dermatophytes in the culture medium depends on the amount and the selected site of suspicious sample material.

4. SPECIMEN COLLECTION AND PREPARATION

- Clean the chosen sampling area with 70 % alcohol to reduce a potential bacterial and/or saprophytic contamination.
- Remove deep scraping from skin and remove hair with root, dandruffs and crusts (with scalpel and clipper) from the edge of the lesions.

5. TEST PROCEDURE

- Spread the obtained material evenly and generously on the agar surface with help of a new sterile scalpel. Leave a small seam in the edge of the tilted agar, so that approaching moulds can be recognized as contaminants.
- Press the spread material well onto the agar surface so that the contact of the hyphae and/or spores to the agar can be ensured. The tighter the contact, the faster the colour changes in the presence of dermatophytes (within 1–3 days) and their growth thereafter.
- Close the lid lightly but not completely and incubate at 25–32°C (77–90°F). Higher temperatures (30–32°C/86–90°F) slow down the growth of mould cultures.
- Check the inoculated sample tubes daily up to 21 days for **colour change** and **colony growth**.
- To avoid false positive test interpretation, it is important to recognize the growth of unwanted saprophyte colonies (mould growth).

6. READING OF THE TEST RESULT

- Colour change:** from orange to red, first after 2–3 days, on average 3–7 days, first indication for dermatophyte growth.
- Colony growth** on average after 5–10 days. Normally, whitish and partially stained colonies appear. In this case, they are yellow-brown at the agar side of the colony (turn sample tube!).

Only a red discoloration around a whitish, softly yellowish to light orange, woolly-fluffy colony refers to the growth of a pathogen dermatophyte.

Mould colonies can be recognized on normally coloured (black-grey-green-brown) growth on the inoculation side of the colonies. Should a colour change of the agar happen, it appears significantly delayed with moulds.

- Macroscopic-visual evaluation** of the colonies:
The identification can be made due to colony size, form and colour as well as their surface structure and texture of the margin.
- Microscopic evaluation** (100-400× magnification):
Take a sample with an inoculation loop from different colony areas and/or at different times. In “cotton-woolly-like” areas rather hyphae can be found, whereas in powdery-plaster-like areas rather spores can be found. The identification results via following criteria: manifestation, width and septation of the hyphae, uniformity of the mycelium, forming of spiral hyphae and/or chlamydospores, formation of micro- (esp. *Trichophyton*) and/or macroconidia (esp. *Microsporum*).

7. PRECAUTIONS FOR USERS

- Make sure that the test-kit components can be correlated exactly to the particular patient.
- Use a new sample tube with tilted agar for each sample.
- The sample material must be seen as potentially infectious and disposed of accordingly, together with the used test-kit components.

8. TEST PRINCIPLE

MYKODERMOASSAY DTM contains a specific selective nutrient agar optimised for the diagnostics of dermatophytes. The dermatophyte medium (DTM) contains specific dermatophyte nutrients, colour indicators and growth inhibiting substances against bacteria and saprophytes (esp. moulds).

The Dermatophyte Test Medium (DTM) belongs to the selective media especially supporting the growth of dermatophytes. Furthermore, it contains phenol red as colour indicator. Dermatophytes particularly use proteins in the first days of growth. The substances resulting thereby cause the colour change of the media after 2–3 days from orange to red. This colour change is an early hint for the growth of dermatophytes in the culture medium.

This enables the veterinarian immediately a fast and targeted initiation or continuation of a specific therapy.

9. INFORMATION FOR THE INTERPRETATION

- The interpretation of the test result should always be based on anamnestic and clinical data as well as the therapy and prophylaxis possibilities.

Culture medium	Colour indicator	Interpretation
DTM (Dermatophyte test medium) Specific dermatophyte medium	Colour change from day 2–3 on orange → red	Growth of dermatophytes → optical differentiation of dermatophyte species

- Pictures of growing dermatophytes on DTM agar can be found on the website http://www.megacor.at/product/mykodermoassay_dtm.html
- In general, in case of obvious mould growth (without/with colour change) the incubation should be stopped because of the rapid and extensive growth of these fungi. Start again with a new sample.
- In general, growth and colour change indicate dermatophytes, exceptionally some moulds (e.g. *Scopulariopsis*, *Chrysosporium* species, “blackness fungi”) can do so, too.
- Treatment should be executed consistently until success of therapy (minimum 6–8 weeks).
- The success of therapy should be controlled with a new culture with MYKODERMOASSAY DTM. In case of a new positive test result, continue with treatment!
- Only two negative culture results at an interval of 4 weeks ensure a therapy success.